

Schizont Transfection using Amaxa Nucleofector – Treeck Lab (June 2026)

Preparation of DNA

The total amount of DNA used per transfection is a minimum of 10 ug of each plasmid (but ideally around 25 ug) in a maximum volume of 10 ul TE. If transfecting a Cas9-sgRNA plasmid with a repair template, make sure to linearise the repair template and use a molar ratio of 1:6 of gRNA plasmid to repair template. Use this website to calculate molarity of your DNA (<https://www.promega.co.uk/resources/tools/biomath/>)

- Precipitate DNA in 2 volumes of 100% EtOH and 0.1 volume 3M Sodium Acetate. You can also precipitate and use the DNA straight away or incubate overnight, but overnight incubation is not necessary.
- Pellet the DNA by spinning at max speed (~20,000 rcf) for 2-5 minutes at room temperature and wash with 700-1000 ul of 70% EtOH twice
- After final EtOH wash, add 700 ul of 70% EtOH, and invert the tube a few times to make sure that the entire internal surface is decontaminated with the EtOH.
- Transfer the tube to a biosafety cabinet and remove as much of the EtOH as possible, then dry the pellet for 15 minutes (residual EtOH will evaporate; make sure to dry completely – DNA will look almost transparent when fully dry)
- Resuspend the DNA in 10 µl TE (sterile) and proceed with the transfection or store in -20°C

Preparation of erythrocytes for invasion after transfection

- Mix 300 ul of packed erythrocytes and 1.5 ml Complete Media in small tissue culture flask (T25) for each sample
- Preheat to 37°C while preparing parasites

Preparation of parasites

*Note 1: Generally, each transfection requires ≥25 ul schizonts. However, using as little as 10 ul schizonts has been successful.

*Note 2: In our experience quantity is more important than quality (i.e., synchronicity).

- In the cycle before, purify schizonts on Percoll cushion and allow the schizonts to invade fresh erythrocytes while shaking. If possible, perform a second synchronisation ~4h later with Sorbitol.
- On the day of transfection, purify schizonts on Percoll cushion once most begin to rupture, then place in RPMI such that there are 25 µl schizonts/500 µl RPMI
- Transfer 500 µl aliquots of resuspended parasites to microcentrifuge tubes
- Pellet parasites briefly in microcentrifuge
- Try to keep parasites warm until you are ready to transfect, since they will slow down and not egress if left at room temperature for too long.

Transfection

- Label AMAXA cuvettes beforehand and pre/open AMAXA pastettes in the biosafety cabinet, since these can be a bit fiddly to open.

- Add 100 ul of supplemented P3 to your DNA that was previously resuspended in 10 ul TE (P3 buffer is supplied with the Lonza transfection kit; beware that this goes off 3 months after supplementation)
- Remove supernatant from pelleted parasites (from above).
- Resuspend parasites in DNA/P3 transfection solution mixture
- Transfer to the appropriately labelled AMAXA cuvette
- Transfect in AMAXA device using setting FP158 for P3 Primary cells
- After transfection, immediately remove cuvette and transfer the parasites with the special AMAXA pastette to a flask with erythrocytes and medium – it is crucial that this is done as quickly as possible.
- After all transfections are done, incubate the flasks with the blood and transfected parasites in the shaking incubator at 37 °C for at least thirty minutes
- Add 8 ml RPMI and place in 37 °C incubator.

Selection and transfection aftercare

The following day, replace medium with medium containing the appropriate drug

Selectable Marker Genes	Drug	Working Concentration	How long for parasites to die fully? (Time kept on drug)	Selection procedure
hDHFR	WR99210 OR Trimethoprim (TMP)	WR: 8x IC ₅₀ (= 0.59 nM) OR TMP: 8x IC ₅₀ (= 10 uM)	4-6 days	Can be added immediately or the day after transfection
BSDr	Blasticidin S deaminase	8x IC ₅₀ (=5 µg/ml)	7-10 days	Add the day after transfection
NeoR	Neomycin (G418)	225 µg/ml	7-10 days	Add the day after transfection
PKG _{T618Q}	ML10	50 nM	7-10 days	Can be added immediately or the day after transfection
yFCU	Ancoti (5-fluoroytosine)	1 mM	2 weeks	Negative selection

Feed transfections every 2-3 days. Make sure to add ~50 ul of fresh blood every week to prevent catastrophic lysis.

Start smearing cultures ~7-10 days post transfection and continue to monitor regularly from this point onwards.